

Portable Fall Alert for Cane, Walker or Wheelchair

This alarm attaches to your walker or your cane. If the walker or cane tips past a 30 degree angle the alarm sounds. This can also be used as a PANIC alert system. Not permanent, easily move it from walker to wheelchair.



ONLY \$48.30 us + Shipping/Handling
*Prices are subject to Change Without Notice
ITEM # ats4525

How do I use this as a FALL ALERT?

Easy, just attach it to the leg of your walker or the post of your cane. If your walker or can falls past a 30° angle, the alarm will go off within 15 seconds. The 15 second delay allows you to pick up the walker or can if you accidentally dropped them. (You can turn off the 15 second delay if you want an instant alert to sound)



How can I use this as a PANIC ALERT?

This alert system has 2 BIG RED buttons that can be used to alert those around you if you need assistance.

One BIG RED button is on the side of the alarm and another BIG RED button that plugs into the fall alarm.

The BIG RED button can be attached to the handle of your cane or walker to be used as a PANIC ALERT.



Just push and hold the BIG RED button and the alarm will keep sounding until you release it. It is as loud as a CAR ALARM - That should get some attention !



ORDER FORM

Assistive Technology Services
"Helping You Live Independently"

8023 Franklin Road
 Murfreesboro, TN, 37128
 (615) 562-0043

Billing Address Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing []	Ship To Address Name _____ Address _____ City _____ State _____ Phone Number _____
---	---	---

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	Ats4525	Portable Fall Alert for Cane, Walker or Wheelchair Shipped to US	\$48.30	\$6.89	
	Ats4525	Portable Fall Alert for Cane, Walker or Wheelchair Shipped to Canada	\$48.30	\$26.89	
	Ats4525	Portable Fall Alert for Cane, Walker or Wheelchair Shipped to Other	\$48.30	\$46.89	
				Subtotal	
		Only Applicable If Shipping To Tennessee		TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice		Shipping	
				Total	

If Paying By Check Please Make Checks Payable To: Assistive Technology Services

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [] MC [] AmEx [] Disc []				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"