

# Steady Spoon Leveling Spoon



## Eat with Confidence

Selection of stabilizing and leveling handles and attachments are designed to help people with hand tremor or limited hand and arm mobility.

**Only \$249.99** us + Shipping/Handling

\*Prices are Subject to Change without Notice

Item# atsss and atsls



## STEADY

### RECHARGEABLE BATTERY

Battery lasts for at least one hour of continuous use

(approximately 3 meals) on a single charge.

Electronic stabilizing handle and a selection of attachments that include a soup spoon, everyday spoon, fork, and spork. Steady is designed to help people with hand tremor, which may be related to Parkinson's disease or essential tremor, eat more easily

It's smart

The stabilizing handle contains sensors that detect hand motion and a small onboard computer that distinguishes unwanted tremor from the intended movement of the hand. To stabilize the utensil, the computer directs two motors in the handle to move the utensil attachment in the opposite direction of any detected tremor.

It's easy to use

Steady includes two parts: a stabilizing handle and a utensil attachment. Connect the attachment to turn it on automatically, enjoy your meal, then simply use the dock to charge.



Connect Attachment



Start Eating



Recharge

## Level

### RECHARGEABLE BATTERY

Battery lasts for at least one hour of continuous use

(approximately 3 meals) on a single charge.

Assistive device for limited hand and arm mobility.

Attached utensil stays level as your hand twists, bends, and moves. Advanced sensor and motor-based self leveling technology.

It's smart

The leveling handle contains sensors that detect changes from the intended movement of the hand in 3 dimensions. To level the utensil, the computer directs two motors in the handle to bend the flexible gray joint, keeping the attachment at the right angle.

Additional utensil attachments can be used with the stabilizing handle including a soup spoon attachment, an everyday spoon attachment, a fork attachment, and a spork attachment.





# ORDER FORM

**Assistive Technology Services**  
*"Helping You Live Independently"*

8023 Franklin Road  
 Murfreesboro, TN, 37128  
 (615) 562-0043  
 Fax (831) 308-9671

<b>Billing Address</b> Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing [ ]	<b>Ship To Address</b> Name _____ Address _____ City _____ State _____ Phone Number _____
-------------------------------------------------------------------------------------------------------	---------------------------------------------	-------------------------------------------------------------------------------------------------------

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING TO US	SHIPPING TO CANADA	SHIPPING TO OTHER	UNIT TOTAL
	ATSSS	Steady Spoon	\$249.99	\$12.89	\$32.89	\$42.89	
	ATSLs	Level Spoon	\$249.99	\$12.89	\$32.89	\$42.89	
	ATSSSXSP	Steady Spoon Spork Attachment	\$49.99	\$8.99	\$18.99	\$28.99	
	ATSSSXES	Steady Spoon Extra Spoon Attachment	\$49.99	\$8.99	\$18.99	\$28.99	
	ATSSSXF	Steady Spoon Fork Attachment	\$49.99	\$8.99	\$18.99	\$28.99	
	ATSSSXSS	Steady Spoon Soup Spoon Attachment	\$49.99	\$8.99	\$18.99	\$28.99	
	ATSLsXSS	Level Spoon Extra Spoon Attachment	\$49.99	\$8.99	\$18.99	\$28.99	
	ATSLsXF	Level Spoon Fork Attachment	\$49.99	\$8.99	\$18.99	\$28.99	
						<b>Subtotal</b>	
		Only Applicable If Shipping To Tennessee				TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice					
						<b>Total</b>	

If Paying By Check Please Make Checks Payable To: Assistive Technology Services

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [ ] MC [ ] AmEx [ ] Disc [ ]				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"