

Talking Wrist Blood Pressure Monitor



This talking, wrist-mounted blood pressure monitor is ideal for those who need to keep track of their blood pressure for medical reasons, and anyone who wants to measure their pulse rate.

With one touch of a button, the monitor auto-inflates for quick, easy measurement, with automatic pressure release valve system and audible data announcements.

ONLY \$57.99 US + Shipping/Handling

*Prices are Subject to Change Without Notice

ITEM # ats4413tw

Measures systolic BP, diastolic BP, mean arterial pressure, and pulse rate.

Memory feature with 2 user-profiles allows you to capture 60 readings with date and time stamp, making it easy to keep track of your results and progress.

Wrist band tucks neatly behind easy-to-read LCD, so the unit can go wherever you go. No more BP-tracking hassle!

General Specifications:

Cuff Type Wrist
Operation Mode Automatic
Talking feature
2 user capability
60 reading memory per user
Date and time stamp
LCD screen size 2" x 1.75"
Measures Systolic, Diastolic, Mean Arterial

Pressure and Pulse Rate

Detailed Specifications:

Number of Users 2 users
Memory 60 memories
Cuff Circumference 5.3" - 7.7"
Measurement Method Oscillo metric method
Measurement Range Pressure: 4.0 to 37.3 kPa (30 to 280 mmHg)
Pulse: 40 to 200 pulse/min
Measurement Accuracy Pressure: +/- 3 mmHg (+/- 0.4 kPa)
Pulse: +/- 5%
Power Source 1 AAA battery
Operating Environment 10°C to 40°C (50°F to 104°F)
Storage Environment -20°C to 60°C (-4°F to 140°F)
Dimensions 3.5" x 3.4" x 1.2"
Weight 6.3 oz.



ORDER FORM

Assistive Technology Services
"Helping You Live Independently"

8023 Franklin Road
 Murfreesboro, TN, 37128
 (615) 562-0043
 Fax (831) 308-9671

Billing Address Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing []	Ship To Address Name _____ Address _____ City _____ State _____ Phone Number _____
---	---	---

CUSTOMER EMAIL ADDRESS		PAYMENT TERMS		DATE
		Due on receipt		

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	Ats4413tw	Talking Wrist Blood Pressure Monitor Shipped to US	\$57.99	\$9.89	
	Ats4413tw	Talking Wrist Blood Pressure Monitor Shipped to Canada	\$57.99	\$29.89	
	Ats4413tw	Talking Wrist Blood Pressure Monitor Shipped to Other	\$57.99	\$49.89	
				Subtotal	
		Only Applicable If Shipping To Tennessee		TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice		Shipping	
				Total	

If Paying By Check Please Make Checks Payable To: Assistive Technology Services

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [] MC [] AmEx [] Disc []				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"